

HARZ Labs Dental IBT

Material Technical Data Sheet (TDS)

Version 1.4 / EN 10 August 2023



SECTION 1: DESCRIPTION AND APPLICATION

Dental IBT is a flexible biocompatible material intended for orthodontic use in order to print indirect bonding trays.

SECTION 2: MATERIAL PROPERTIES

2.1 Characteristics of liquid

Tested property	Standard/Method	Result
Color	-	Clear
Odor	-	Weak
Density	ASTM D1298	1.1 ± 0.1 g/cm ³
Viscosity (25 °C)	ASTM D2393	550 ± 150 mPa*s

2.2 Mechanical properties

Tested property	Standard/Method	Result
Ultimate Tensile Strength	ASTM D638	20.0 ± 5.0 MPa
Tensile Modulus	ASTM D638	250.0 ± 35.0 MPa
Elongation at Break	ASTM D638	120 ± 20 %
Hardness	ASTM D2240	75 ± 5 Shore D

2.3 Special parameters

2.0 openial parameters		
Tested property	Standard/Method	Result
Solubility (24h)	ASTM D3132	≤ 0.01 %
Sorption (24h)	ASTM D570	≤ 0.47 %
Free Monomer Extraction	ISO 10993	Passed
Cytotoxicity	ISO 10993	Passed
Irritation	ISO 10993	Passed
Sensitization	ISO 10993	Passed

The information above is believed to be accurate and represents the best information currently available to us. All test specimens were printed, cleaned, and post-processed per instructions provided by HARZ Labs company. Results provided here are representative of these processes and may vary if these established protocols are not followed. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall HARZ Labs LLC (OOO «XAPLĮ Лабс») be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if HARZ Labs LLC (OOO «XAPLĮ Лабс») has been advised of the possibility of such damages.